







Bone Cutting Forceps




Ruskin-Liston Forceps
D 537 400 6" Straight




Kleinert-Kutz Forceps
Delicate Blades
D 537 450 6" Straight
D 537 470 6" Curved




Ruskin-Liston Forceps
Straight
D 537 280 7-1/2"
D 537 300 7-1/2" Delicate Blades



Bone Cutting Forceps



Stille-Liston Forceps
D 535 560
10-1/2" Straight




Stille-Liston Forceps
D 535 600
10-1/2" Angled on Flat




Stille-Horsley Forceps
D 535 640
10-1/2" Angled

Rongeurs




Micro Friedman Rongeur
1.5x5mm Bite

- D 534 300 5-1/2" Straight
- D 534 320 5-1/2" 30° Angled
- D 534 340 5-1/2" 90° Angled




Friedman Rongeur

- 4x11mm Bite
- D 534 520
- 5-1/2" Curved Delicate



Mini-Friedman Rongeur


- 2x8mm Bite, Delicate
- D 534 530 6" Straight
- D 534 540 6" Curved



Dean (Blumenthal) Rongeur

- 4x11mm Bite
- D 534 560 6" 30° Curve

Rongeurs



1/1
D 534 480




1/1
D 534 500

Lempert Rongeur

3x7mm Bite, Delicate

D 534 480 6" Straight

D 534 500 6" Curved



1/1
D 534 440




1/1
D 534 460

Lempert Rongeur

3x7mm Bite, Delicate

D 534 440 7-1/2" Straight

D 534 460 7-1/2" Curved



1/1
D 536 480



1/1
D 536 500



Stille-Luer Rongeur

10x18mm Bite

D 536 480 8-1/2" Straight

D 536 500 8-1/2" Curved


Rongeurs and Trephines



Leksell Rongeur

9", 15° Curve

D 536 820	3x16mm bite	D 536 780	6x16mm bite
D 536 760	4x16mm bite	D 536 790	7x16mm bite
D 536 770	5x16mm bite	D 536 800	8x16mm bite




*Designed for bone biopsy.
Also facilitates removal of
broken bone screws.*



Michele Trephine

D 544 290	7" Small, 3mm inside diameter, 5mm outside diameter
D 544 300	7" Medium, 4mm inside diameter, 6mm outside diameter
D 544 310	7" Large, 6mm inside diameter, 8mm outside diameter



1/1








P.I. Plug Cutter

D 544 420

4-1/2" with Obturator

3mm inside diameter, 5mm outside diameter



- These cutters are especially useful for arthrodesis of small and large joints by plug cutter graft technique.
- Plug cutters are designed to be used in pairs.
- Smaller diameter instrument will cut the plug opening and the larger one will harvest a bone graft that will fit into the opening.



P.I. Plug Cutter

6"

D 760 200 5mm inside diameter, 7mm outside diameter

D 760 220 7mm inside diameter, 9mm outside diameter

D 760 240 10mm inside diameter, 12mm outside diameter

D 760 260 12mm inside diameter, 14mm outside diameter

D 760 280 20mm inside diameter, 22mm outside diameter

D 760 300 22mm inside diameter, 24mm outside diameter