

Mallets




Collin Mallet

D 528 802
8" Solid Stainless Steel
Weight: 7 oz. (200 grams)
Head weight: 5 oz. (150 grams)
Head diameter: 30mm



Lucae Bone Mallet

D 528 500
8" Brass Head
Stainless Steel Handle
Weight: 9 oz. (260 grams)
Head Weight: 5 oz. (150 grams)
Head Diameter: 25mm



Combination Mallet

D 528 760
7-1/2"
1 Nylon Cap and 1 Solid Stainless Cap
Weight: 9 oz. (260 grams)
Head Weight: 5 oz. (150 grams)
Head Diameter: 26mm

Nylon Caps Autoclavable to 250°F (121°C). Replacement caps available, contact your sales representative for details.



Nylon Mallet

D 528 520
7-1/2" Lightweight
2 Replacable Nylon Caps
Anodized Aluminum Handle
Weight: 7 oz. (200 grams)
Head Weight: 4 oz. (120 grams)
Head Diameter: 25mm

