

# Mallets




### Collin Mallet

D 528 802  
8" Solid Stainless Steel  
Weight: 7 oz. (200 grams)  
Head weight: 5 oz. (150 grams)  
Head diameter: 30mm



### Lucae Bone Mallet

D 528 500  
8" Brass Head  
Stainless Steel Handle  
Weight: 9 oz. (260 grams)  
Head Weight: 5 oz. (150 grams)  
Head Diameter: 25mm



### Combination Mallet

D 528 760  
7-1/2"  
1 Nylon Cap and 1 Solid Stainless Cap  
Weight: 9 oz. (260 grams)  
Head Weight: 5 oz. (150 grams)  
Head Diameter: 26mm

*Nylon Caps Autoclavable to 250°F (121°C). Replacement caps available, contact your sales representative for details.*



### Nylon Mallet

D 528 520  
7-1/2" Lightweight  
2 Replacable Nylon Caps  
Anodized Aluminum Handle  
Weight: 7 oz. (200 grams)  
Head Weight: 4 oz. (120 grams)  
Head Diameter: 25mm