



Wire and Pin Cutters




1/1

Cuts flush to the bone like no other wire cutter, a big advantage in buried wire procedures.




P.I. Flush Wire Cutter

D 528 210
 7" Double Action, Front and Side Cutting
 Carbide Inserts
 Maximum capacity .062" (1.6mm)




1/1

1/1


The smooth double action combined with the versatile jaw design make this cutter a favorite with many surgeons.

Double Action Wire Cutter

D 528 130
 7" Angled
 Side Cutting with Carbide Inserts
 Maximum capacity .062" (1.6mm)



1/1



Cannulation allows you to slide the cutter over the inserted wire. The piece cut off will be held safely inside the channel.

Cannulated Pin Cutter

D 528 370
 7-1/2" End Cutting
 Maximum capacity .062" (1.6mm)