



## Drill Bit JC

D 543 190  
127mm length, 1.0mm diameter  
D 543 200  
127mm length, 1.5mm diameter  
D 543 220  
127mm length, 2.0mm diameter  
D 543 240  
127mm length, 2.4mm diameter  
D 543 260  
127mm length, 2.7mm diameter

D 543 270  
127mm length, 3.2mm diameter  
D 543 280  
127mm length, 3.5mm diameter  
D 543 300  
127mm length, 4.0mm diameter  
D 543 320  
127mm length, 4.7mm diameter  
D 543 330  
127mm length, 6.0mm diameter

*Designed to fit any  
Jacobs chuck drill.*