


Osteotomes



The unusually long bevel of this osteotome allows the cut bone to gradually bend as it is cut.



P.I. Beveled Osteotome


7" Straight

- D 761 300 10mm
- D 761 360 20mm
- D 761 400 25mm

7" Curved

- D 761 600 5mm
- D 761 650 10mm

Osteotomes



Hoke Osteotome

5-1/4" Straight

- D 526 870 1/8" (3mm)
- D 526 880 3/16" (4mm)
- D 526 890 1/4" (6mm)
- D 526 900 5/16" (8mm)
- D 526 910 3/8" (10mm)
- D 526 920 1/2" (13mm)
- D 526 930 5/8" (16mm)
- D 526 940 3/4" (19mm)

5-1/4" Curved

- D 528 900 1/8" (3mm)
- D 528 910 3/16" (4mm)
- D 528 920 1/4" (6mm)
- D 528 930 5/16" (8mm)
- D 528 940 3/8" (10mm)
- D 528 960 1/2" (13mm)
- D 528 980 5/8" (16mm)
- D 528 990 3/4" (19mm)