







# K-Wire Driver and Gratloch™ Wire Bender



- Drive, bend and pull K-Wire, all in one tool
- One handed operation
- Powerless, ideal for office or hospital use
- Contains depthscale
- No set-up or maintenance
- Fully autoclavable
- For wire from .028" to .062"

## K-Wire Driver

D 527 300  
5-1/2" German Stainless



*The Gratloch™ Wire Bender stabilizes and bends K-Wire simultaneously. Wire can be bent close to the bone single-handedly without fracture or osteotomy stability.*



## P.I. Gratloch™ Wire Bender

D 545 600  
7-1/2" Hardened, Wear-resistant Stainless Steel  
Maximum capacity .062" (1.6mm)  
U.S. Patent #5,314,431