


# Super-Cut Bandage Scissors


**Super-Cut Scissors** represent an advance in scissor technology. Created from the finest surgical steels with master craftsmanship, the upper blade has a razor-sharp edge which gives the scissors a smooth and effective cutting action. However, they are not recommended for fiberglass or plaster splints. Additional benefits:

- Less hand fatigue for the surgeon and less trauma for the patient.
- Retain their sharpness longer than ordinary scissors.
- Ring handles are black-coated for easy identification and are fully autoclavable.




## Super-Cut Hi-Level Bandage Scissors

D 413 920 5-1/2" Angled Shanks  
Delicate offset shanks and sharp Super-Cut edge provide better leverage. The fine knob slides easily under finger and toe bandages.



## Super-Cut Lister Bandage Scissors with Large Finger Ring

D 414 040 8"  
All the benefits of our Super-Cut Lister scissors plus an oversized finger ring, which helps relieve pressure on the third metacarpal.



## Super-Cut Lister Bandage Scissors

D 414 000 5-1/2"  
D 414 020 7-1/4"


# Suture and Bandage Scissors



### Shortbent Stitch Scissors


D 412 700  
3-1/2" Curved, Delicate

*Our Stitch Scissors feature delicate hooked blades for cutting suture.*



### Spencer Stitch Scissors

D 412 580  
3-1/2" Straight, Delicate  
D 412 820  
4-3/4" Straight, Delicate



1/1  
D 378 940

### Lister Bandage Scissors

D 378 940  
7-1/4"  
D 378 980  
8" One Large Finger Ring