





Self-Retaining Retractors



These Retractors are designed to provide excellent visibility. Each is self-retaining and automatically expands the wound site.


Low Profile Small Wound Retractor

D 586 160 2-1/2" x 3"
 Curved Smooth Blades
 10mm depth x 12mm height
 Maximum opening: 30mm
 Center blade: 10mm x 16mm



Low Profile 2 Prong Retractor

D 586 180 2-1/2" x 3"
 2 x 2 Blunt Prongs
 8mm depth x 17mm height
 Maximum opening: 30mm
 Center blade: 10mm x 16mm



Low Profile 3 Prong Retractor

D 586 170 3" x 3-3/4"
 3 x 3 Blunt Prongs
 11mm depth x 25mm height
 Maximum opening: 35mm
 Center blade: 10mm x 16mm