

Mallets




Collin Mallet

D 528 802
8" Solid Stainless Steel
Weight: 7 oz. (200 grams)
Head weight: 5 oz. (150 grams)
Head diameter: 30mm



Lucae Bone Mallet

D 528 500
8" Brass Head
Stainless Steel Handle
Weight: 9 oz. (260 grams)
Head Weight: 5 oz. (150 grams)
Head Diameter: 25mm



Combination Mallet

D 528 760
7-1/2"
1 Nylon Cap and 1 Solid Stainless Cap
Weight: 9 oz. (260 grams)
Head Weight: 5 oz. (150 grams)
Head Diameter: 26mm

Nylon Caps Autoclavable to 250°F (121°C). Replacement caps available, contact your sales representative for details.



Nylon Mallet

D 528 520
7-1/2" Lightweight
2 Replacable Nylon Caps
Anodized Aluminum Handle
Weight: 7 oz. (200 grams)
Head Weight: 4 oz. (120 grams)
Head Diameter: 25mm


Mallets and Bone Tamps




Short Ortho Mallet
 D 559 530
 7" Lightweight
 Solid Stainless Steel, Short Handle
 Weight: 1 lb. 2 oz. (500 grams)
 Head weight: 8 oz. (240 grams)
 Head diameter: 30mm



Ortho Mallet
 D 528 801
 9-1/2" Aluminum Handle
 Weight: 1 lb. 2 oz. (510 grams)
 Head diameter: 30mm/40mm




1/1



Non-slip grid design allows for tapping bone grafts or bone parts into place with minimal trauma to the bone.



P.I. Bone Impactor
 D 761 700
 6-1/2"
 Anodized Aluminum Handle
 Replaceable Nylon Caps
 Weight: 1 lb. 2 oz. (500 grams)
 Head diameter: 20mm



Bone Tamp

Cross Serrated Tip

- D 530 830 6", 2mm diameter
- D 530 840 6", 3mm diameter
- D 530 850 6", 4mm diameter
- D 530 860 6", 5mm diameter
- D 530 870 6", 6mm diameter
- D 530 880 6", 8mm diameter
- D 530 890 6", 10mm diameter