





*Delicate needle nose jaws specifically designed to handle cerclage wire. Jaw tapers to 1mm.*

**Needle Nose Pliers**  
D 544 320  
5-1/2" Delicate





**Needle Nose Pliers**  
D 544 350  
5-1/2", 2mm tip



*Combination plier and cutter. Cuts all sizes (0.2mm to 1.5mm) of cerclage wire. Jaw tapers to 3.5mm.*

**Flat Nose Pliers**  
D 544 330  
6-1/2", 3.5mm tip



**Flat Nose Pliers**  
D 544 360  
5-1/2" Broad, 5mm tip