

# Mallets and Bone Tamps




## Short Ortho Mallet

D 559 530  
 7" Lightweight  
 Solid Stainless Steel, Short Handle  
 Weight: 1 lb. 2 oz. (500 grams)  
 Head weight: 8 oz. (240 grams)  
 Head diameter: 30mm




## Ortho Mallet

D 528 801  
 9-1/2" Aluminum Handle  
 Weight: 1 lb. 2 oz. (510 grams)  
 Head diameter: 30mm/40mm



1/1




*Non-slip grid design allows for tapping bone grafts or bone parts into place with minimal trauma to the bone.*



## P.I. Bone Impactor

D 761 700  
 6-1/2"  
 Anodized Aluminum Handle  
 Replaceable Nylon Caps  
 Weight: 1 lb. 2 oz. (500 grams)  
 Head diameter: 20mm



## Bone Tamp

Cross Serrated Tip

D 530 830 6", 2mm diameter  
 D 530 840 6", 3mm diameter  
 D 530 850 6", 4mm diameter  
 D 530 860 6", 5mm diameter  
 D 530 870 6", 6mm diameter  
 D 530 880 6", 8mm diameter  
 D 530 890 6", 10mm diameter