



Nail Splitters

- Our nail splitters have the true feel of superb craftsmanship. The cutting edge is expertly ground and specially heat treated and the grooved handles provide a sure grip.
- Nail splitters should not be used on thick or mycotic nails, use nail nippers for these applications.



This short jaw design provides increased leverage for easier cutting of thick nails, causing less strain on the hand and splitter.





P.I. Short Jaw Nail Splitter

Stainless

D 761 200 5"

D 761 220 6"

D R11 330 Replacement Springs (2)



Nail Splitter

Stainless

D 755 100 4"

D 755 120 4-1/2"

D 755 160 5" Medium Jaw

D 755 170 5" Narrow Jaw

D 755 300 6"

D R11 330 Replacement Springs (2)

D R11 440 Replacement Springs (2)
(for D 755 100 only)