


# Pliers and Plate Bending Press




1/1

1/1

*Combination pliers and wire cutter has all the benefits of our Parallel Pliers plus a cutting jaw for K-Wires up to .062" (1.6mm).*

## Parallel Plier with Cutter

D 547 950  
7", 10mm tip  
Maximum opening 1/2"  
Maximum capacity .062" (1.6mm)



1/1

## Mini Plate Bending Pliers

D 544 340  
5" for 1.5mm and 2.0mm plates



## Plate Bending Press

D 541 190  
12"