




# Skin Hooks and Retractors



**Frazier Hook**  
 D 543 480  
 5" Sharp  
 D 543 510  
 5" Blunt



**Kleinert-Kutz Hook**  
 D 582 500  
 6" 3mm diameter  
 D 582 540  
 6" 5mm diameter  
 D 582 580  
 6" 7mm diameter




**Freer Skin Hook**  
 D 582 300  
 6" 2 Sharp Prongs,  
 2.5mm wide  
 D 582 320  
 6" 2 Blunt Prongs,  
 2.5mm wide



**Joseph Hook**  
 D 642 780  
 6-1/4" 1 Sharp Prong  
 D 642 810  
 6-1/4" 2 Sharp Prongs,  
 2mm wide



**Volkman Finger Retractor**  
 D 583 880  
 4-1/2" 2 Prongs, Sharp  
 D 583 980  
 4-1/2" 2 Prongs, Blunt



**Rigid Rake Retractor**  
 2 Prongs  
 D 583 480 6-1/2" Sharp  
 D 583 560 6-1/2" Blunt