



# Nail Splitters

- Our nail splitters have the true feel of superb craftsmanship. The cutting edge is expertly ground and specially heat treated and the grooved handles provide a sure grip.
- Nail splitters should not be used on thick or mycotic nails, use nail nippers for these applications.



*This short jaw design provides increased leverage for easier cutting of thick nails, causing less strain on the hand and splitter.*





### P.I. Short Jaw Nail Splitter

Stainless

D 761 200 5"

D 761 220 6"

D R11 330 Replacement Springs (2)



### Nail Splitter

Stainless

D 755 100 4"

D 755 120 4-1/2"

D 755 160 5" Medium Jaw

D 755 170 5" Narrow Jaw

D 755 300 6"

D R11 330 Replacement Springs (2)

D R11 440 Replacement Springs (2)  
(for D 755 100 only)