

# Mallets and Bone Tamps




## Short Ortho Mallet

D 559 530  
7" Lightweight  
Solid Stainless Steel, Short Handle  
Weight: 1 lb. 2 oz. (500 grams)  
Head weight: 8 oz. (240 grams)  
Head diameter: 30mm



## Ortho Mallet

D 528 801  
9-1/2" Aluminum Handle  
Weight: 1 lb. 2 oz. (510 grams)  
Head diameter: 30mm/40mm




*Non-slip grid design allows for tapping bone grafts or bone parts into place with minimal trauma to the bone.*



## P.I. Bone Impactor

D 761 700  
6-1/2"  
Anodized Aluminum Handle  
Replaceable Nylon Caps  
Weight: 1 lb. 2 oz. (500 grams)  
Head diameter: 20mm



## Bone Tamp

Cross Serrated Tip  
D 530 830 6", 2mm diameter  
D 530 840 6", 3mm diameter  
D 530 850 6", 4mm diameter  
D 530 860 6", 5mm diameter  
D 530 870 6", 6mm diameter  
D 530 880 6", 8mm diameter  
D 530 890 6", 10mm diameter