


# Super-Cut Scissors



## Super-Cut Iris Scissors

Sharp Points


- D 416 400 3-1/2" Straight, Delicate
- D 416 420 3-1/2" Curved, Delicate
- D 416 440 4-1/2" Straight
- D 416 500 4-1/2" Curved

## Super-Cut Stevens Scissors

Blunt Points

- D 416 740 4-1/2" Straight
- D 416 800 4-1/2" Curved
- D 416 840 5-1/4" Curved

*Fine sharp blades for cutting delicate tissue during plastic, hand, and minor vascular surgery.*



## Super-Cut Operating Scissors

Sharp/Blunt Points

- D 414 100 5-1/2" Straight
- D 414 400 5-1/2" Curved

## Super-Cut Jamison Scissors

- D 416 860 6-1/4" Curved, Blunt Points
- Very Delicate blunt scissors with long handles for fine tissue dissection