


Super-Cut Bandage Scissors

Super-Cut Scissors represent an advance in scissor technology. Created from the finest surgical steels with master craftsmanship, the upper blade has a razor-sharp edge which gives the scissors a smooth and effective cutting action. However, they are not recommended for fiberglass or plaster splints. Additional benefits:

- Less hand fatigue for the surgeon and less trauma for the patient.
- Retain their sharpness longer than ordinary scissors.
- Ring handles are black-coated for easy identification and are fully autoclavable.



Super-Cut Hi-Level Bandage Scissors

D 413 920 5-1/2" Angled Shanks
Delicate offset shanks and sharp Super-Cut edge provide better leverage. The fine knob slides easily under finger and toe bandages.


Super-Cut Lister Bandage Scissors with Large Finger Ring

D 414 040 8"
All the benefits of our Super-Cut Lister scissors plus an oversized finger ring, which helps relieve pressure on the third metacarpal.

Super-Cut Lister Bandage Scissors

D 414 000 5-1/2"
D 414 020 7-1/4"

Super-Cut Scissors



Super-Cut Iris Scissors

Sharp Points


- D 416 400 3-1/2" Straight, Delicate
- D 416 420 3-1/2" Curved, Delicate
- D 416 440 4-1/2" Straight
- D 416 500 4-1/2" Curved

Super-Cut Stevens Scissors

Blunt Points

- D 416 740 4-1/2" Straight
- D 416 800 4-1/2" Curved
- D 416 840 5-1/4" Curved

Fine sharp blades for cutting delicate tissue during plastic, hand, and minor vascular surgery.



Super-Cut Operating Scissors


Sharp/Blunt Points

- D 414 100 5-1/2" Straight
- D 414 400 5-1/2" Curved

Super-Cut Jamison Scissors

- D 416 860 6-1/4" Curved, Blunt Points
- Very Delicate blunt scissors with long handles for fine tissue dissection

Super-Cut Scissors




1/1
D 415 540




Super-Cut Baby Metzenbaum Scissors

Blunt Points

- D 415 540 4-1/2" Straight
- D 415 600 4-1/2" Curved




1/1
D 415 720



Super-Cut Metzenbaum Scissors

Blunt Points


- D 415 700 6" Straight
- D 415 720 6" Curved
- D 415 820 7" Straight
- D 415 840 7" Curved
- D 415 860 7" Curved, Delicate
- D 415 960 8" Curved
- D 416 080 9" Curved



Super-Cut Mayo Scissors

Beveled Blades, Blunt Points

- D 414 700 5-1/2" Straight
- D 414 820 6-3/4" Straight
- D 415 060 9" Straight
- D 414 760 5-1/2" Curved
- D 414 880 6-3/4" Curved
- D 415 120 9" Curved



1/1
D 414 700



1/1
D 414 820



1/1
D 415 060



1/1
D 415 120