


A new instrument available exclusively from Delta Surgical Instruments!

The Modified Inge Retractor with K-Wire Guide



- Accommodates up to .062 K-Wire.
- Can be used as a distractor and/or retractor.
- Ratchet graduated up to 35mm to accurately assess how wide the lateral portion of the bone graft needs to be.
- Useful in many types of procedures to accurately assess amount of bone graft needed, such as fusions and nonunions.
- Narrow tips allow for easy entrance into small incisions.
- Serrated tips for non-traumatic tissue abrasion.



 **Delta**
Surgical Instruments

To place an order call Toll Free 866-390-2226 • www.deltasurgicalinstruments.com