

A new instrument available exclusively from Delta Surgical Instruments!

# The Fibular Sesamoid Elevator



## The precise instrument for an effective lateral release in Hallux Valgus Surgery

- Narrow and curved shaft  
Allows selective and precise release, minimizing dissection and exposure
- Concave cutting edge  
Captures capsular layer within the tip – cuts only *within* the tip
- Contoured to match the shape of the metatarsal head and the fibular sesamoid  
Allows quick but safe access to specific tissues to be released



**T**he *Fibular Sesamoid Elevator*\* is the right instrument; designed specifically to facilitate quick and precise release of the fibular sesamoid portion of a lateral release. It allows the surgeon to precisely access and release soft tissues limiting motion or correction without excessive exposure.

 **Delta**  
Surgical Instruments



Developed in conjunction with Thomas F. Smith, DPM  
and *The Podiatry Institute*.

\*Patent Pending

To place an order call Toll Free 866-390-2226 • [www.deltasurgicalinstruments.com](http://www.deltasurgicalinstruments.com)